# GIMA

PROFESSIONAL MEDICAL PRODUCTS

## **OCCHIALI INGRANDENTI BINOCULAR LOUPES LOUPES BINOCULAIRES LUPAS BINOCULARES LUPAS BINOCULARES ΜΕΓΕΘΥΝΤΙΚΕΣ ΔΙΟΠΤΡΕΣ LOUPES**

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All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

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Zumax Medical Co., Ltd.  
5 Zhiying Street, Suzhou New District,  
Suzhou, Jiangsu, China 215129 - Made in China

MedNet EC-REP GmbH Borkstrasse 10,  
48163 Münster, Germany

Gima S.p.A. Via Marconi, 1 - 20060 Gessate (MI) Italy  
gima@gimaitaly.com - export@gimaitaly.com  
[www.gimaitaly.com](http://www.gimaitaly.com)

**PRODUCT CARE**

- Do not expose your product to liquid or moisture or to humidity.
- Do not touch the surface of lens with hand or hard objects.
- Never clean the surface of lens with any hard and dirty materials.
- Do not look into the sun or into a powerful light source with this loupe.

## ***Purpose and features of the product***

### **Purpose**

SL series Binocular Loupes is designed for all possible medical needs. They are widely accepted by Stomatology, ENT, Gynecology, Cardiology, Brain surgery and Plastic surgery etc..

### **Features**

- Multi-coated, high grade glass lenses produce crisp, clear images.
- The Loupes which employ eyepiece of very large view field diameter. The depth of focus is deep and object can be seen stereoscopically with actual distance feeling.
- The angle of the loupes can be adjusted between 0°~40° with ease and fatigue is minimized even after prolonged use.
- 6mm up and down displacement. Pupil distance can be adjusted individually which could be suited for all kind of faces.
- Adapts to any Head band and Headlight of Zumax
- Smart design, convenient for mobile diagnosis.

## ***Operating Instruction***

First step, please fasten the strap to the frame of glasses.

Put on the loupe and check that the pupil distance is symmetrically set.

Now look at a well it three-dimensional object. Adjust one eyepiece after the other until you see a round image with no shadows and a focused object in each. Then screw in both set screws.

The position of the optics is adjusted over a wide range by means of the hinge mount. The ease of movement of the hinges can be adjusted by means of the hinge screws.

## ***Maintenance of product***

### **Maintenance of lens**

- The instrument not in use should be covered and put into the packing box.
- Do not disassemble the optical lens by unskilled or unqualified technicians.

### **Cleaning**

- The outer surface of the instrument may be cleaned with wet cloth, the remaining stains can be cleaned off with the mixture of 50% alcohol and 50% distilled water. Do not wipe with any corrosive detergent lest that the surface should be damaged.
- Cleaning the lenses: Mirror-cleaning paper or a drop of liquid solvent (50% alcohol and 50% aether) may be used to clean the lens, blow it carefully. If there is dust stained on the lens, blow them with a blowball or brush them with a dust pen.

**Caution: Be careful to prevent the solvent into the instrument**

Specifications

| Magnification | Working distance (mm) | Field of view dia. (mm) |
|---------------|-----------------------|-------------------------|
| 2.5x          | 340/420               | 80/90                   |
| 3.5x          | 340/420               | 55/60                   |

Environment requirements

|           |                         |                    |
|-----------|-------------------------|--------------------|
| Operation | Environment Temperature | +10°C – +40°C      |
|           | Relative humidity       | 30% – 75%          |
|           | Atmospheric pressure    | 700 hPa – 1060 hPa |

|          |                   |
|----------|-------------------|
| Shipping | Common conveyance |
|----------|-------------------|

|         |                         |                    |
|---------|-------------------------|--------------------|
| Storage | Environment Temperature | +10°C – +40°C      |
|         | Relative humidity       | 30% – 75%          |
|         | Atmospheric pressure    | 700 hPa – 1060 hPa |

|                                                                                   |                                                 |                                                                                   |                              |                                                                                   |                                                        |
|-----------------------------------------------------------------------------------|-------------------------------------------------|-----------------------------------------------------------------------------------|------------------------------|-----------------------------------------------------------------------------------|--------------------------------------------------------|
|   | Product code                                    |  | Lot number                   |  | Medical Device                                         |
|   | Keep in a cool, dry place                       |  | Keep away from sunlight      |  | Manufacturer                                           |
|   | Caution: read instructions (warnings) carefully |  | Consult instructions for use |  | Medical Device compliant with Regulation (EU) 2017/745 |
|   | Date of manufacture                             |  | Imported by                  |  | Unique device identifier                               |
|  | Model number                                    |                                                                                   |                              |                                                                                   |                                                        |

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.