GIMA

PROFESSIONAL MEDICAL PRODUCTS

BRACCIALE PER INFUSIONE
INFUSION CUFF
BRASSARD DE PERFUSION A PRESSION
INFUSIONSMANSCHETTE
BRAZAL PARA INFUSIÓN
MANGUEIRA DE INFUSÃO
ΕΓΧΥΜΑΤΙΣΜΕΝΟ ΜΑΝΣΑΚΙ
تسريب الكفة

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All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

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Maclin Power Ltd
20 Wenlock Road London N1 7GU
England, United Kingdom

Infusion Cuff Set

This device is applicable to be used as the ancillary of infusion in clinic, and it is recommended to be used together with the infusion bag. Comparing with the traditional infusion method with pressure provided by gravity, this device can avoid the contamination of the ambient air as no ambient air will contact with the liquid medicine. In the meantime, it makes the paramedic more convenient as the paramedic don't need to hold the infusion set at given level accompanying the patient, especially when the patients go for eating, outdoor activity or toilet.

The device consists of infusion cuff, inflation bulb, manometer and knot. Please connect these parts before use.

The instruction of using this instrument

This device can provide necessary pressure for the infusion bag so as to ensure the liquid medicine flows from the infusion bag to the patients' body.

- Put the infusion bag into the infusion cuff, and the infusion bag will be just below the transparent net. Put the infusion cuff in safe area to avoid extrusion or falling down, fix the infusion bag by the pressure of infusion cuff;
- Pump the inflation bulb to inflate air into the inflation cuff. Observing the manometer carefully to make sure that the pressure reaches the average of the patient's systolic pressure and diastolic pressure;
- Connect the disposable infusion set with the infusion bag, and the liquid medicine will flow into the whole disposable infusion set. Make sure that there is not bubble in it, and close the infusion switch. Then it is ready for other ordinary infusion operation;
- During the infusion, the pressure will decrease gradually. Please pump the inflation bulb to increase the pressure. If the infusion speed is too low or too high, please modify the pressure in the infusion cuff by pumping the bulb to inflate or releasing deflation valve to deflate so as to get expected infusion speed;

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

REF	IT Codice prodotto GB Product code FR Code produit DE Erzeugniscode ES Código producto PT Código produto GR Κωδικός προϊόντος SA كود المنتج	MD	IT Dispositivo medico GB Medical Device FR Dispositif médical DE Medizinprodukt ES Producto sanitario PT Dispositivo médico GR Ιατροτεχνολογικό προϊόν SA جهاز طبي
LOT	IT Numero di lotto GB Lot number FR Numéro de lot DE Chargennummer ES Número de lote PT Número de lote GR Αριθμός παρτίδας SA رقم الدفعة		IT Leggere le istruzioni per l'uso GB Consult instructions for use FR Consulter les instructions d'utilisation DE Gebrauchsanweisung beachten ES Consultar las instrucciones de uso PT Consulte as instruções de uso GR Διαβάστε προσεκτικά τις οδηγίες χρήσης SA اقرأ بدقة وحرص تعليمات الاستخدام
	IT Conservare in luogo fresco ed asciutto GB Keep in a cool, dry place FR À conserver dans un endroit frais et sec DE An einem kühlen und trockenen Ort lagern ES Conservar en un lugar fresco y seco PT Armazenar em local fresco e seco GR Διατηρείται σε δροσερό και στεγνό περιβάλλον SA يحفظ في مكان بارد وجاف		IT Conservare al riparo dalla luce solare GB Keep away from sunlight FR À conserver à l'abri de la lumière du soleil DE Vor Sonneneinstrahlung geschützt lagern ES Conservar al amparo de la luz solar PT Guardar ao abrigo da luz solar GR Κρατήστε το μακριά από ηλιακή ακτινοβολία SA يحفظ بعيداً عن ضوء الشمس
	IT Fabbricante GB Manufacturer FR Fabricant DE Hersteller ES Fabricante PT Fabricante GR Παραγωγός SA الشركة المصنعة		IT Data di fabbricazione GB Date of manufacture FR Date de fabrication DE Herstellungsdatum ES Fecha de fabricación PT Data de fabrico GR Ημερομηνία παραγωγής SA تاريخ التصنيع
	IT Dispositivo medico conforme al regolamento (UE) 2017/745 GB Medical Device compliant with Regulation (EU) 2017/745 FR Dispositif médical conforme au règlement (UE) 2017/745 DE Medizinprodukt im Sinne der Verordnung (EU) 2017/745 ES Producto sanitario conforme con el reglamento (UE) 2017/745 PT Dispositivo médico em conformidade com a regulamentação (UE) 2017/745 GR Ιατρική συσκευή σύμφωνα με την ΚΑΝΟΝΙΣΜΟΣ (ΕΕ) 2017/745 (UE) 2017/745 SA جهاز طبي يتوافق مع التوجيه		IT Contenuto o presenza di lattice naturale GB Contains or presence of natural rubber latex FR Contient ou présence de latex de caoutchouc naturel DE Enthält oder Anwesenheit von Naturlatex ES Contiene o presencia de látex de caucho natural PT Contém ou presença de látex de borracha natural GR Περιέχει ή περιέχει λατέξ φυσικού καουτσούκ SA يحتوي على أو وجود مطاط طبيعي من المطاط
UDI	IT - Identificatore unico dispositivo GB - Unique device identifier FR - Identifiant unique de l'appareil DE - Eindeutige Gerätekennung ES - Identificador de dispositivo único PT - Identificador exclusivo do dispositivo GR - Μοναδικό αναγνωριστικό συσκευής SA - معرف الجهاز الفريد	CH REP	IT - Rappresentante autorizzato in Svizzera GB - Authorized Representative in Switzerland FR - Représentant autorisé en Suisse DE - Bevollmächtigter Vertreter in der Schweiz PT - Representante autorizado na Suíça GR - Εξουσιοδοτημένος Αντιπρόσωπος στην Ελλάδα SA - الممثل المعتمد في سويسرا
UK REP	IT - Rappresentante autorizzato nel Regno Unito GB - Authorized Representative in the UK FR - Représentant autorisé au Royaume-Uni DE - Bevollmächtigter Vertreter in Großbritannien ES - Representante autorizado en el Reino Unido PT - Representante autorizado no Reino Unido GR - Εξουσιοδοτημένος Αντιπρόσωπος στο Ηνωμένο Βασίλειο SA - الممثل المعتمد في المملكة المتحدة		

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Gima S.p.A.